AE

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User Instructions (GB)

Cotton Aragó

Cotton band
Sterile and non-resorbable surgical ligation

Description:
Cotton Aragó is a sterile, twisted and non-resorbable surgical ligation comprised of multi-filament twisted cotton 30/2 prepared by spinning and twisting of the cotton fibres from Languo Gossypil absorbens. Available in white with no optical dyes or whiteners. Meets the requirements for hydrophilic cotton according to European Pharmacopoeia and American USP.

Indications:
Cotton Aragó is used as ligation for the majority of tissues.

Effect:
Cotton Aragó is used as a ligation and doesn't remain in the body.

Contraindications:
Unknown.

Warnings:
The user of Cotton Aragó must be familiar with surgical techniques and the use of non-resorbable ligations. Do not resterilise/reuse. Do not use if the packaging is open. Reusing a single-use product presents a potential infection risk for the patient. Contamination of the product may cause an injury, a disease or death for the patient. Reprocessing may compromise the structural integrity of the product and may cause failure of said product. The absence of labelling and instructions at the point of use may result in the product being used after its expiry date and losing its traceability. Modification of the original packaging may compromise the sterility and functionality of the product.

Precautions:
The handling of this or any other ligation requires care. Avoid damaging the ligation through crushing or folding it by using tweezers. The surgeon should consider the type of knotting used in the surgical techniques. When handling the ligation, appropriate precautions should be taken to avoid surgical instruments such as tweezers damaging the ligation material due to crushing, folding or abrasion. The user must be familiar with surgical ligation techniques.

Side effects:
Not applicable.

Presentations:
This ligation is supplied in strands precut to different lengths and with a variable number of strands. Consult the catalogue for more detailed information.

Storage conditions:
Store at room temperature, avoid sustained exposure to high temperatures and protect against humidity. Do not sure beyond the expiry date.

Explanation of symbols on the packaging:

LOT	Batch No
REF	Catalogue reference
	Single use only
	Do not re-sterilise
	Expiry yyyy-mm- year/month
	See instructions
STERILE R	Sterilised by Gamma Rays
	Keep out of sunlight
	Do not use if the unit packaging is open or deteriorated
	Manufacturer
CE 0318	CE Marking + identification number of the notified body. This product meets the essential requirements of the Council Directive 93/42/EEC concerning medical devices

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