(AR)

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Trade name of the product (EN)

PDO Aragó

Polydioxanone Monofilament
Surgical, sterile and absorbable suture

Manufacturer's name and registered office
Laboratorio Aragó S.L.
Salvador Mundi 11
08017 Barcelona (Spain).

Description/General information
PDO Aragó is a sterile, synthetic, monofilament, absorbable surgical suture composed of a homopolymer of p-dioxanone. The empirical molecular formula is $(C_4H_6O_2)_n$.
For better visibility in the surgical field, the thread is dyed with D&C Violet No. 2 dye (Colour index number CI 60725), although it can be supplied undyed.
PDO Aragó suture is available in a wide variety of gauges and lengths and are packaged with or without a stainless steel atraumatic needle of various types and sizes, as described in the PRESENTATION section. Check what the device includes on the label on the outside.
PDO Aragó complies with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile, Synthetic, Absorbable, Monofilament Sutures and the United States Pharmacopeia (USP) for Synthetic Absorbable Sutures. The European Pharmacopoeia recognises the equivalence of the american and Ph. Eur. units, which is detailed on the labelling.
PDO Aragó sutures are designed only to be used by healthcare professionals trained in surgical suturing techniques. They are intended to be used in patients for any of the indications included in the device's indications for use. It is not indicated for elderly people, those who are frail or have slow-healing wounds.
The expected clinical benefits of PDO Aragó are:
- Holding the edges of a wound or injury in apposition (together) to aid the natural healing process, so reducing dehiscence, infection and pain.
- Use as a ligature to occlude vessels (arteries, veins or vascularised tissues) and thus mechanically prevent bleeding.
The summary of safety and clinical performance can be found at the following link (after activation): https://ec.europa.eu/tools/eudamed

Intended purpose
PDO Aragó is an absorbable suture composed of a homopolymer of p-dioxanone. Used in surgical procedures for approximation of the edges of a soft tissue wound, providing temporary support for the wound until it reaches a degree of healing that allows it to withstand normal tension, or as a ligature for vessel occlusion to prevent haemorrhage. Used by healthcare professionals in patients of all age ranges. Not suitable for cardiovascular or neurological surgery.

Indications
PDO Aragó is indicated for applications where, in accordance with surgical practice, absorbable materials are used for sutures. It is suitable for the approximation of soft tissue and ligatures, but not for cardiovascular or neurological tissue.

Contraindications
Since PDO Aragó is absorbable, it must not be used on wounds that depend on prolonged reinforcement with sutures. It is not indicated for elderly people, those who are frail or have slow-healing wounds.
The safety and efficacy of PDO Aragó suture in cardiovascular or neurological tissues have not been established.

Warnings
Reusing this product (or any of its parts) carries a risk of product degradation, which may result in its failure or in cross-contamination, which in turn may lead to infection or transmission of bloodborne pathogens to patients and healthcare professionals.
Mild reactions may appear initially in tissue around the suture material. Due to its absorbable nature, this suture may act as a foreign body temporarily. Prolonged contact of a suture with saline solutions, such as those found in the urinary or biliary tract, may cause stones to form. Like all foreign bodies, sutures can aggravate an existing infection.
Any processing may compromise the structural integrity of the product or may cause it to fail. Absence of labelling and instructions at point of use may result in the device being used after its expiry date and traceability being lost. Modification of the original packaging may compromise the sterility and functionality of the device.
Care should be taken when using absorbable sutures in tissues with poor blood supply, as suture extrusion and delayed absorption may occur.
Adverse reactions associated with the use of this suture include tissue separation or wound dehiscence leading to poor healing, or failure to provide adequate wound support upon closure of the sites where expansion, stretching or distension occur. In tissues with poor blood supply, suture extrusion may occur.
This also occurs in elderly, malnourished or debilitated patients, or in patients with other conditions that may delay wound healing.
Breakage of the suture may cause hemorrage. Suture material that cuts or tears the tissue on tensioning the suture can causes lesions in soft tissue. Improper knotting or damage to the suture during use may result in prolonged surgery or failure of the treatment and further surgery.
Needle breakage may result in prolonged surgery, require additional surgery or become residual foreign bodies. Accidental punctures with contaminated surgical needles can lead to the transmission of blood-borne pathogens.

Precautions
Before using PDO Aragó sutures for wound closure, healthcare professionals should be familiar with the surgical procedures and techniques in which sutures are used, as the risk of wounds reopening (dehiscence), leading to impaired healing, depends on the site of use and the suture material used.
Proper surgical practice should be followed when treating contaminated or infected wounds.
Since this is an absorbable suture, healthcare professionals should consider using a supplementary non-absorbable suture when suturing areas that may be subjected to expansion, tension or distension, or that require additional support.

As with any suture material, good knot fixation requires the application of the standard surgical technique of flat and square ties with additional loops depending on the surgical circumstances and the surgeon's experience. Improper knotting or damage to the suture during use may result in prolonged surgery or failure of the treatment and further surgery.
Cutaneous sutures that need to remain in place for more than 7 days may cause local irritation and should therefore be cut or removed as directed by a healthcare professional.
Subcutaneous sutures should be placed as superficially as possible to minimise the erythema and induration normally associated with the absorption process.
This suture may be unsuitable for elderly, malnourished or weak patients, or those with conditions that delay wound healing.
When handling sutures, care should be taken to ensure that surgical instruments must be taken to prevent surgical instruments such as forceps and needle holders from damaging the suture material by crushing, folding or scraping and/or deforming the needle.
The needle should be held in an area one third (1/3) to one half (1/2) of the distance from the distal end to the needle tip. If it is held close to the tip area, it will reduce its penetration capacity and may also cause the needle to fracture. If the needle is held at its distal end, it may cause it to bend or break.
If needles are reshaped manually they may lose strength and be less resistant to bending or breaking.
Healthcare professionals must be familiar with surgical suturing techniques and must handle needles carefully in order to avoid injuries caused by accidental punctures, that can transmit blood-borne pathogens from contaminated needles.
In case of malfunction of the device prior to use, e.g. a bent, broken or detached needle, or a damaged suture, the device should be discarded and a new one obtained to start the procedure. If the device fails to work properly during use, the healthcare professional will decide whether or not to continue with the product and to stop or to start a new procedure.
The effect that chemical and/or radiological treatments may have on suture functionality must be taken into account.
Healthcare professionals must inform the patient of the adverse reactions, undesirable side effects and risks associated with the device and the procedure, and recommend the patient contact a healthcare professional in the event of abnormal post-operative recovery.
Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority in the country.

Application/Instructions for use
Sutures should be selected and implanted according to the patient's condition, surgical experience, surgical technique and the nature of the wound.
Dispose of used needles in a sharps container. Dispose of unintentionally opened, partially used or used devices and their packaging in accordance with your institution's policies and procedures for bio-hazardous materials and waste.
The accredited organisation is responsible for the safe disposal of the device. Any unused devices or waste material to be disposed of in accordance with local requirements for potentially contaminated materials.
Radiological or chemical treatments may only be carried out once the wounds have healed. At this point, the suture is encapsulated or hydrolysed and no longer performs its function.

Effect / Performance / Mechanism of action
PDO Aragó causes a minimal inflammatory tissue reaction with gradual encapsulation due to connective tissue. Gradual loss of resistance to breakage and eventual absorption of PDO Aragó occurs through hydrolysis, producing glycolate, which is then absorbed and metabolised by the body. Absorption starts with loss of tensile strength, followed by loss of mass.
Implantation studies in rats show the following profile: After 2 weeks, the suture retains approximately 70-80% of the initial tensile strength, depending on the USP number. At 4 weeks remaining tensile strength is 60-70%, with absorption completed in approximately 180 - 220 days.

Sterility
PDO suture is sterilised using ethylene oxide. Do not re-sterilise. Do not use if packaging is open or damaged.

Storage and handling conditions
Store at room temperature, avoid sustained exposure to high temperatures. Protect from humidity and direct sunlight. Do not use after expiry date.

Presentations
This suture can be supplied in sizes USP 7/0 (EP 0,5) to USP 2 (EP 5). It is supplied sterile, in pre-cut strands of different lengths and variable numbers of strands, without needle or threaded with different types of stainless-steel atraumatic needles of various sizes. The suture is available in boxes of 12 or 36 units. Consult the catalogue for more detailed information.

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