(AE)

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User Instructions (GB)

PDO H Aragó

Several Polydioxanone filaments intertwined
Sterile, absorbable surgical suture

Description:
PDO H Aragó is a sterile, synthetic, absorbable surgical suture composed of several monofilaments of Polydioxanone (PDO) intertwined. The empirical formula is $(C_4H_6O_3)_n$. For better visibility, the thread has been coloured with D&C Violet No 2 (C.I.60725) colouring. It meets the requirements of the European Pharmacopoeia for Sterile Synthetic Absorbable Monofilament Sutures and of the American Pharmacopoeia USP for Synthetic Absorbable Sutures, except in the diameter.

Indications:
PDO H Aragó is indicated for use in surgical procedures where absorbable suture materials are used. It is indicated for use in soft tissue approximation and ligation. It is not indicated for use in skin, paediatric and/or adult cardiovascular tissue, ophthalmic surgery, microsurgery or neurological tissue.

Effect:
PDO H Aragó elicits a minimal initial reaction in tissue with gradual encapsulation by connective tissue. Gradual loss of resistance to breakage and eventual absorption of PDO H Aragó occurs through hydrolysis. The polymer transforms to 2-hydroxyethoxy acetic acid which is subsequently absorbed and metabolised by the body. Absorption starts with a loss of tensile strength followed by mass loss. Implantation studies in rats show the following profile:
After 2 weeks, the suture retains approximately 70-80 % of the initial tensile strength, depending on the USP number. At 4 weeks remaining tensile strength is 60 - 70% with absorption completed in about 180 - 220 days.

Contraindications:
Because PDO H Aragó is absorbable it must not be used on wounds that require prolonged reinforcement with suture.

Warnings:
Users of PDO H Aragó should be familiar with surgical techniques using absorbable sutures because the risk of wound reopening (dehiscence) may vary with the site of application and the suture material used. The use of PDO H Aragó may not be appropriate for elderly or weak patients or those whose wounds heal slowly. Because it is an absorbable suture, the surgeon must consider the use of supplemental nonabsorbable suture for places that may undergo expansion, stretching or distension or which may require additional support. Do not resterilise/reuse. Do not use if the packaging is open. Reusing a single-use product presents a potential infection risk for the patient. Contaminated product may cause injury, disease or death of the patient. Reprocessing may compromise the structural integrity of the product and cause it to fail. Absence of labelling and instructions at point of use may result in the product being used after its expiry date and traceability being lost. Changes to the original packaging may compromise the sterility and functionality of the product.

Precautions:
Skin sutures that must remain in place for more than 7 days may cause localised irritation and should be snipped off or removed. In tissues with deficient blood supply the thread may be rejected due to slow absorption. This or any other suture must be handled carefully. Avoid crushing or crimping damage to the suture due to the use of forceps or needle holders. If needles are reshaped manually they may lose strength and be less resistant to bending or breaking. The surgeon should consider the type of knotting used in the surgical techniques. When handling sutures, appropriate precautions must be taken to prevent surgical instruments such as forceps and needle holders from damaging the suture material by crushing, crimping or abrasion and/or from deforming the needle. Users of these sutures must be familiar with surgical suturing techniques and must handle needles carefully to avoid accidental needle sticks. Dispose of used needles in a dedicated container. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the distal end to the point. Grasping the needle towards the point will reduce its penetration capacity and may also cause the needle to fracture. Grasping the needle at the distal end may cause it to twist or fracture.

Side effects:
Initially slight reactions may appear in tissue around the suture material. Like any foreign body, PDO H Aragó may aggravate a prior infection. As with any foreign body, prolonged contact of any suture with salt solutions such as those found in the urinary and biliary tracts may result in calculus formation.

Presentations:
This suture is supplied with the following composition PDO H Aragó in precut threads of different lengths, without needle or threaded with different types of atraumatic needle. See the catalogue for more detailed information.

Storage conditions:
Store at room temperature, avoid sustained exposure to high temperatures and protect against humidity. Do not use beyond the expiry date.

Explanation of symbols on the packaging:

LOT	Batch No.
REF	Catalogue reference
	Single use only.
	Do not re-sterilise.
	Expiry yyyy-mm year/month.
	See instructions.
STERILE EO	Sterilised with ethylene oxide.
	Keep out of sunlight.
	Do not use if the unit packaging is open or damaged.
	Manufacturer
CE 0318	EC Marking + identification number of the notified body. This product meets the essential requirements of Council Directive 93/42/EEC concerning medical devices.

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