AE

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User Instructions (GB)

# Polypropylene Aragó

## Surgical, sterile and non-absorbable suture

**Description:**
Polypropylene Aragó is a sterile, monofilament and non-absorbable surgical suture based on Polypropylene. For better visibility, the thread has been dyed blue with Solvent Blue 104 (C.I. 61568) dye. It complies with the requirements of Sterile Non-Absorbable Fibres according to the European and American (USP) Pharmacopoeias.

**Indications:**
Polypropylene Aragó is indicated for applications where, in accordance with surgical practice, non-absorbable materials are used for sutures. It is suitable for the approximation of soft tissue and for ligatures, including use in cardiovascular and neurological tissue, ophthalmology, nerve repairs and aponeurosis, plastic surgery and abdominal walls.

**Effect:**
Polypropylene Aragó causes a minimal acute reaction of the tissue on which a gradual encapsulation occurs due to the connective tissue. It is resistant to traction, slides easily through tissue, and has high chemical and biological inertness. It is non-absorbable and remains in the body indefinitely.

**Contraindications:**
None known.

**Warnings:**
Polypropylene Aragó users must be familiar with surgical techniques and the use of non-absorbable sutures, because the risk of the injury reopening (dehiscence) depends on where it is used and the suture material employed. Do not re-sterilise/re-use. Do not use if the packaging is open. Re-using a single-use product creates a potential risk of infection for the patient. Product contamination may cause injury or disease, or lead to the death of the patient. Re-processing could jeopardise the structural integrity of the product or cause it to fail. Using the product without the label and instructions will result in using the product beyond its expiry date and the loss of its traceability. Altering the original packaging may jeopardise the sterility and functionality of the product.

**Precautions:**
All sutures must be handled with care. Avoid damaging the suture by crushing or folding it through the use of forceps or needle holders. If the needles are remodelled manually, they may lose their strength and become less resistant to bending and breaking. The surgeon must take into consideration the type of knot to be used in a surgical technique. When handling sutures, due precautions must be taken to prevent surgical instruments such as forceps and needle holders from damaging the suture material by crushing, folding or scraping it, and/or deforming the needle. The user must be familiar with suture techniques and must handle the needles with care so as to not be pricked by them. Discard used needles into a special container. Hold the needle between one-third (1/3) and a half (1/2) of the distance from the distal extremity to the tip of the needle. If it is held close to the tip, the capacity for penetration will decrease and this could also cause a fracture of the needle. If the needle is held by the distal extremity, the needle may bend or fracture.

**Side effects:**
Mild reactions in the tissue surrounding the suture material may appear at the beginning. As with any foreign body, Polypropylene Aragó may exacerbate a prior infection. As with any foreign body, prolonged contact of this or any other suture with saline solutions such as those found in urinary and biliary conduits may lead to the formation of calculus.

**Presentations:**
This suture can be supplied in USP sizes 10/0 (EP 0.2) up to USP 2 (EP 5) in pre-cut threads of different lengths and with a variable number of threads, without a needle or swaged to different types of atraumatic needles. See the catalogue for more detailed information.

**Storage conditions:**
Store at room temperature, avoid prolonged exposures to high temperatures and protect from dampness. Do not use after the expiry date.

**Explanation of symbols on the packaging:**

| | |
|---|---|
| LOT | Batch No |
| REF | Catalogue reference |
| | Single use only |
| | Do not re-sterilise |
| | Expiry yyyy-mm- year/month |
| | See instructions |
| STERILE EO | Sterilised with ethylene oxide |
| | Keep out of sunlight |
| | Do not use if the unit packaging is open or deteriorated |
| | Manufacturer |

CE 0318 CE Marking + identification number of the notified body. This product meets the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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