AR

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Effet/Performance/Mécanisme d'action
La suture PTFE Aragó suture causes a minimal inflammatory tissue reaction on which a gradual encapsulation occurs due to the connective tissue.

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Instructions for use (EN)

PTFE Aragó

Polytetrafluoroethylene monofilament
Sterile non-absorbable surgical suture

Name of the manufacturer and the address of its registered place of business
Laboratorio Aragó S.L.
Salvador Mundi 11
08017 Barcelona-España (Spain).

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Intended purpose
PTFE Aragó is a non-absorbable monofilament suture of polytetrafluoroethylene. Used in surgical procedures for approximation of the edges of a soft tissue wound, providing temporary support for the wound until it reaches a degree of healing that allows it to withstand normal tension, or as a ligature for vessel occlusion to prevent haemorrhage. Used by healthcare professionals in patients of all age ranges. Not to be used for use in cardiovascular, ophthalmic or neurological surgery.

Indications
PTFE Aragó is indicated for application where non-absorbable suture materials are used for surgical procedures. It is suitable for soft tissue approximation and/or ligatures, including use in dental surgery, as well as general surgery, but not for use in cardiovascular surgery, ophthalmology, microsurgery or in direct contact with the central nervous system.

Contraindications
The safety and efficacy of PTFE Aragó sutures in cardiovascular, ophthalmic or neurological tissues has not been established.

Warnings
Re-use of this device (or any of its parts) carries a risk of cross-contamination, which may lead to infections or transmission of blood-borne pathogens to patients and healthcare professionals.
At first, mild reactions may appear in the tissue around the suture material. Like any foreign body, PTFE Aragó may aggravate a previous state of infection. Prolonged contact with this or any other suture with saline solutions such as those from urinary and biliary tracts can lead to the formation of calculi.
Reprocessing may jeopardise the structural integrity of the product or may cause it to fail. Using the product without the label and instructions will result in using the product beyond its expiry date and the loss of its traceability. Altering the original packaging may jeopardise the sterility and functionality of the product.
Adverse reactions associated with the use of this suture include tissue separation or wound dehiscence, which leads to poor healing or failure to provide adequate wound support upon closure of the sites where expansion, stretching or distension occur. In tissues with poor blood supply, suture extrusion may occur.
Suture breakage can result in haemorrhage. Suture material cutting or tearing through the tissue at the time of suture tensioning may cause soft tissue injury. Improper knotting or damage caused to the suture during use may result in additional surgery, or the failure of treatment and further surgery.
Broken needles may result in extended or additional surgeries or residual foreign bodies. Accidental punctures with contaminated surgical needles can lead to the transmission of blood-borne pathogens.

Precautions
Before using PTFE Aragó suture for wound closure, healthcare professionals should be familiar with the procedures and surgical techniques in which sutures are used, as the risk of tissue separation (dehiscence) may lead to impaired healing, depending on the site and the suture material used.
Proper surgical practice should be followed for the management of contaminated or infected wounds.
Like any suture material, adequate knot fixation requires the standard surgical application of square flat and square ties with additional loops as warranted by surgical circumstances and experience of the surgeon. Improper knotting or damage caused to the suture during use may result in additional surgery, or the failure of treatment and further surgery.
When handling sutures, proper precautions should be taken to prevent surgical instruments, such as forceps and needle holders from damaging the suture material by crushing, folding, or scraping and/or deforming the needle.
The needle should be held in an area between one-third (1/3) and one-half (1/2) of the distance from the distal end to the needle tip. If it is held close to the tip area, the capacity of penetration will decrease and this could also cause a fracture of the needle. If the needle is held at its distal end, the needle may bend or fracture.
If needles are manually re-shaped, they may lose their strength and be less resistant to bending and breaking.
Healthcare professionals should be familiar with surgical suturing techniques and should handle needles with caution to avoid injuries caused by accidental punctures that can transmit blood-borne pathogens from contaminated needles.
Should a product malfunction before use, for example a bent, broken or detached needle or a damaged suture, the product should be discarded and a new one obtained to begin the procedure. Document the event. If a product not working properly during its use, it is at the discretion of the healthcare professional to continue or discontinue usage of the product and how to complete the procedure.
The effect that chemical and/or radiological treatments may have on the functionality of the suture must be taken into account.
Healthcare professionals should inform the patient of adverse reactions, undesirable side effects and risks associated with the product and the procedure and advise the patient to contact a healthcare professional in case of any change in the normal post-operative course.
Any serious incident that has occurred in relation to the product must be reported to the manufacturer and the competent authority in the country.

Application/ Instructions for use
Sutures should be selected and implanted according to the patient's condition, surgical experience, surgical technique and wound characteristics.
Dispose of used needles in a sharps container. Dispose of unintentionally opened, partially used, or used products and their packaging in accordance with your facility's policies and procedures for biohazardous materials and waste.
The accredited organisation is responsible for safe disposal of the product. Any product or waste material must be disposed of in accordance with local requirements for potentially contaminated materials.

**Effect /Performance/M