(AR)

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Description/ General information
Supramid Aragó is a sterile, synthetic, non-absorbable surgical suture. It can have different configurations depending on its diameter:
- USP 4/0 to 6: pseudomonofilament formed by twisted polyamide 6/6 (C12H22N2O2) filament coated with a polyamide 6 (NH-(CH2)5-CO) sheath, giving it a monofilament appearance.
- USP 7/0 to 5/0: Monofilament in polyamide 6 and/or polyamide 6/6.

For better visibility, the thread is coloured black with Hematein/Logewood Black (CI 75290) or Pigment Black 7 (CI 77266) dyes. It can also be supplied colourless.
Supramid Aragó sutures are available in a wide variety of gauges and lengths, and are packaged with or without stainless-steel atraumatic needles of various types and sizes, as described in the PRESENTATION FORMATS section. Check what the device includes on the label on the outside.
Supramid Aragó sutures comply with the requirements of the European Pharmacopoeia (Ph. Eur.) for Sterile Non-Absorbable Synthetic Sutures and the American Pharmacopoeia (USP) for Non-Absorbable Sutures. The European Pharmacopoeia recognises the equivalence of the American and Ph. Eur. units, which is detailed on the labelling.
Supramid Aragó sutures are designed only to be used by healthcare professionals trained in surgical suturing techniques. They are intended to be used in patients for any of the indications included in the device's indications for use.
The expected clinical benefits of Supramid Aragó are:
- Holding the edges of a wound or injury in apposition (together) to aid the natural healing process, so reducing dehiscence, infection and pain.
- Use as a ligature to occlude vessels (arteries, veins or vascularised tissues) and thus mechanically prevent bleeding.

The summary of safety and clinical performance can be found at the following link (after activation): https://ec.europa.eu/tools/eudamed

Intended use
Supramid Aragó is a twisted non-absorbable suture made of polyamide. Used for the approximation of wound edges in soft tissue for natural healing or as a ligature to occlude vessels and prevent bleeding. Used by healthcare professionals in patients of all age ranges. Not suitable for use in cardiovascular, ophthalmic or neurological surgery.

Indications
Supramid Aragó is indicated for application where non-absorbable suture materials are used for surgical procedures. It is suitable for soft tissue approximation and for ligatures, except in cardiovascular, ophthalmic or neurological tissue.

Contraindications
The safety and efficacy of Aragó Supramid sutures in cardiovascular or neurological tissues has not been established.

Warnings
Re-use of this device (or any of its parts) carries a risk of cross-contamination, which may lead to infections or transmission of blood-borne pathogens to patients and healthcare professionals.
Mild reactions may initially appear in tissue around the suture material. Like any foreign body, Supramid Aragó may aggravate a prior state of infection. Prolonged contact of this or any other suture with saline solutions such as those found in the urinary and biliary tracts may result in calculus formation.
Reprocessing may compromise the structural integrity of the product or may cause it to fail. Absence of labelling and instructions at point of use may result in the device being used after its expiry date and traceability being lost. Modification of the original packaging may compromise the sterility and functionality of the device.
Adverse reactions associated with the use of this suture include tissue separation or wound dehiscence leading to poor healing, or failure to provide adequate wound support upon closure of the site where expansion, stretching or distension occur. In tissues with poor blood supply, suture extrusion may occur.
Breakage of the suture may cause haemorrhage. Suture material that cuts or tears the tissue on tensioning of the suture can cause lesions in soft tissue. Improper knotting or damage to the suture during use may result in prolonged surgery or failure of the treatment and further surgery.
Needle breakage may result in prolonged surgery, require additional surgery or become residual foreign bodies. Accidental punctures with contaminated surgical needles can lead to the transmission of blood-borne pathogens

Precautions
Before using Supramid Aragó sutures for wound closure, healthcare professionals should be familiar with the surgical procedures and techniques in which sutures are used, as the risk of wounds reopening (dehiscence), leading to impaired healing, depends on the site of use and the suture material used.

Proper surgical practice should be followed when treating contaminated or infected wounds.
As with any suture material, good knot fixation requires the application of the standard surgical technique of flat and square ties with additional loops depending on the surgical circumstances and the surgeon's experience. Improper knotting or damage to the suture during use may result in prolonged surgery or failure of the treatment and further surgery.
When handling sutures, appropriate precautions must be taken to prevent surgical instruments such as forceps and needle holders from damaging the suture material by crushing, folding or scraping and/or deforming the needle.
The needle should be held in an area one third (1/3) to one half (1/2) of the distance from the distal end to the needle tip. If it is held close to the tip area, it will reduce its penetration capacity and may also cause the needle to fracture. If the needle is held at its distal end, it may cause it to bend or fracture.
If needles are reshaped manually they may lose strength and be less resistant to bending or breaking.
Healthcare professionals must be familiar with surgical suturing techniques and must handle needles carefully in order to avoid injuries caused by accidental punctures, that can transmit blood-borne pathogens from contaminated needles.
In case of malfunction of the device prior to use, e.g. a bent, broken or detached needle, or a damaged suture, the device should be discarded and a new one obtained to start the procedure. If the device fails to work properly during use, the healthcare professional will decide whether or not to continue with the procedure and how to complete the procedure.
The effect that chemical and/or radiological treatments may have on suture functionality must be taken into account.
Healthcare professionals should inform the patient of the adverse reactions, undesirable side effects and risks associated with the device and the procedure, and recommend the patient contact a healthcare professional in the event of abnormal post-operative evolution.
Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority in the country.

Application/Instructions for use
Sutures should be selected and implanted according to the patient's condition, surgical exper