

Hand Tendon Kit

Supramid / Nylon

Supramid Loop



Supramid Loop 4/0

Nylon 6/0 ATB-10

Nylon Skin TB-12

2 X Supramid Loop 4/0

Kessler Technique (1st sachet)

Tang Technique (2nd sachet)

3/8 circle needle round

REF	NEEDLE	LENGTH OF THE THREAD	COLOUR	USP	METRIC
B12	16 mm	30 cm	black	4/0	1,5

3/8 circle needle cutting precision point

REF	NEEDLE	LENGTH OF THE THREAD	COLOUR	USP	METRIC
ATB10	13 mm	22 cm	black	6/0	0,7

3/8 circle needle cutting

REF	NEEDLE	LENGTH OF THE THREAD	COLOUR	USP	METRIC
TB12	16 mm	45 cm	black	4/0	1,5

Hand Tendon Kit 19801

Packaging in boxes of 2 units

Check available variations in thread length and needle with Laboratorio Aragó


Supramid

Hand Tendon Kit

Nylon

Structure	Twisted multifilament coated with a sheath of polyamide that gives it the appearance and behaviour of a monofilament	Monofilament
Composition	Polyamide 6 and Polyamide 6/6	
Colour	Black	
Sterilization	Gamma irradiation	
Absorption	Non absorbable	
Characteristics	<ul style="list-style-type: none"> •Excellent tensile strength •Very good manageability due to its multifilament and monofilament structure •Easy slipping through tissues •Easiness in knot tying •Minimal memory effect 	<ul style="list-style-type: none"> •Excellent tensile strength •Excellent manageability due to its low memory effect •High knot security •Easy slipping through tissues
Indications	Kit specific for the professional of tendon-hand surgery	
Packaging	Boxes of 2 units	

Supramid Loop USP 4/0 HTK Aragó
Polypropylene USP 4/0 Pharmacopoeia



- LOT Batch number
- REF Catalogue reference
- Do not reuse
- Do not resterilise
- Expiry date yyyy-mm
- STERILE R Sterilized by Gamma Rays
- Keep away from sunlight
- Do not use if package is damaged
- See instructions of use
- Manufacturer

Supramid
CE 0318 CE-mark identification number of notified body. Product fulfills requirements of the medical Device Directive 93/42/EEC

Nylon
CE 0318 CE-mark identification number of notified body. Product fulfills requirements of the medical Device Directive 93/42/EEC



Aragó Sutures

Hand Tendon Kit

Supramid / Nylon

We facilitate the functional
recovery of the patient



Laboratorio Aragón

Salvador Mundí, 11
08017 Barcelona
Tel. 93 204 22 12
Fax 93 203 26 09
export@laboratorioarago.com
www.laboratorioarago.com